

# Lightweight Handheld Detachable Compliant Robotic Laryngoscope with Lightweight Intelligent Visual Guidance

Tao Zhang<sup>1</sup>, Sattar Kadir <sup>2</sup>, Hongyu Geng <sup>1</sup>, Huilin Pan <sup>1,3</sup>, An Wang <sup>1</sup>, Jiewen Lai <sup>1,\*</sup> and Hongliang Ren <sup>1,\*</sup>


- <sup>1</sup> Department of Electronics Engineering, The Chinese University of Hong Kong
- <sup>2</sup> Department of Mechanical and Automation Engineering, The Chinese University of Hong Kong
- <sup>3</sup> Department of Artificial Intelligence, Huazhong University of Science and Technology



# 1 Background \_\_\_


### What is tracheal intubation(TI)?

Tracheal intubation is the act of passing a tube into the trachea via the mouth or nose (although access to the trachea can also be created via a direct surgical approach, namely a cricothyrotomy or tracheostomy).



**Current approach** 

### How to insert the tracheal tube using our robot?





**Higher dexterity!** 

Lower difficulty!

Al-based guidance!

## 2 Mechanical design





# 2 Kinematic modeling

#### **Forward kinematics**

$$l_i = nl_i + (-1)^i \cdot \frac{\theta d_i}{2}$$
, given  $i \in \{1, 2\}$ 

$$\theta = (-1)^i \cdot \frac{2(l_i - nl_i)}{d_i},$$

$$\begin{cases} {}_{1}^{0}T = \operatorname{Trans}(P_{10})\operatorname{Rot}(z,\phi)\operatorname{Rot}(y,\theta) \\ {}_{2}^{1}T = \operatorname{Trans}(P_{21}) \\ {}_{2}^{0}T = {}_{1}^{0}T \cdot {}_{2}^{1}T \end{cases}$$

$$\mathcal{R}_z( heta_z) = egin{bmatrix} \cos heta_z & -\sin heta_z & 0 \ \sin heta_z & \cos heta_z & 0 \ 0 & 0 & 1 \end{bmatrix} \; \mathcal{R}_y( heta_y) = egin{bmatrix} \cos heta_y & 0 & \sin heta_y \ 0 & 1 & 0 \ -\sin heta_y & 0 & \cos heta_y \end{bmatrix}$$

#### **Inverse kinematics**

$$\phi = \arctan\left(\frac{P_y}{P_x}\right)$$

$$P_z = \frac{nl_i}{\theta} \sin \theta + l_t \cos \theta$$


$$\lim_{\theta \to 0} \frac{\sin \theta}{\theta}$$


$$\sin \theta = x - \frac{\theta^3}{3!} + \frac{\theta^5}{5!} + o(\theta^7)$$

$$\cos \theta = 1 - \frac{\theta^2}{2!} + \frac{\theta^4}{4!} + o(\theta^6)$$

$$P_z = nl_i \left( 1 - \frac{\theta^2}{3!} + \frac{(\theta^2)^2}{5!} \right) + l_t \left( 1 - \frac{\theta^2}{2!} + \frac{(\theta^2)^2}{4!} \right)$$


### A lightweight solution






# 4 Experiment

#### Verification of the kinematic model:




Experimental set up

### **Evaluation of the glottis detection:**




The phantom



The experiment result indicates that the position error is smaller than 2.34 mm, which is 2.7% over the total length of the flexible segment and the fiducial marker.


Experiment results




Under the blurry condition, our detector could identify the glottis with over 70% confidence coefficient and a dynamic bounding box, with a model size of less than 5 MB and a time cost of around 50 ms.

#### The detection results

# 4 Experiment





Pig carcass experiment set up



In vitro view



In vivo view

## 5 Main contribution & future work \_\_\_\_\_

### The core contributions of this work include the following:

- Lightweight Compliant Robotic Laryngoscope: An all-in-one lightweight detachable compliant robotic laryngoscope (0.65 kg) is developed. The robotic laryngoscope comprises a steerable flexible segment with an endoscopic camera and a three-DoF (bending, axial rotation, and anchoring) flexible segment that can be steerable in one hand.
- **Detachability:** A detachable design is adopted to allow the disassembly and assembly of the steerable flexible segment swiftly to ensure sanitation.
- **Lightweight Intelligent Visual Guidance:** A self-contained, lightweight glottis detection network displays the target glottis at the laryngoscope monitor in real time is proposed.

#### **Future work:**

- Structure optimization of the steerable segment and ergonomic optimization for the system.
- Development of more Al-based surgical process guidance functions



## Thanks for listening!

Best Regards



Scan for more about our lab!